

Occlude a wide range of vessel sizes with coil lengths from 3-20 cm and diameters from 2-20 mm.



Now available in more sizes.

**Nester**<sup>®</sup>  
EMBOLIZATION COIL

**Expand your options with a large number of coil profiles**

Use with a complementary line of embolization products to support your procedure from start to finish.

**Excellent packing coil**

Soft platinum coils pack densely with soft fiber for excellent occlusion.

**Use alone or in combination**

Use in conjunction with other pushable coils, vascular plugs and detachable coil systems.



[www.cookmedical.com](http://www.cookmedical.com)

## NOW AVAILABLE IN A BROADER SELECTION OF MICROCOILS

Rows shaded in blue are new additions to the product line.

	Order Number	Reference Part Number	Recommended Catheter ID/ End Hole Diameter inch	Extended Embolus Length cm	Coiled Embolus Diameter mm	Approximate Number of Loops
<b>.018 Inch Microcoils™</b>						
3cm Extended Embolus	G52731	MWCE-18-3-2-NESTER	.018	3	2	4.8
	G52732	MWCE-18-3-3-NESTER	.018	3	3	3.2
5cm Extended Embolus	G52733	MWCE-18-5-2-NESTER	.018	5	2	8.0
	G52734	MWCE-18-5-3-NESTER	.018	5	3	5.3
	G23520	MWCE-18-5-5-NESTER	.018	5	5	3.2
7cm Extended Embolus	G52735	MWCE-18-7-2-NESTER	.018	7	2	11.1
	G52736	MWCE-18-7-3-NESTER	.018	7	3	7.4
	G52737	MWCE-18-7-4-NESTER	.018	7	4	5.6
	G23521	MWCE-18-7-5-NESTER	.018	7	5	4.5
	G52738	MWCE-18-7-6-NESTER	.018	7	6	3.7
	G52739	MWCE-18-7-8-NESTER	.018	7	8	2.8
	G52741	MWCE-18-7-10-NESTER	.018	7	10	2.2
14cm Extended Embolus	G26987	MWCE-18-14-3-NESTER	.018	14	3	14.9
	G26988	MWCE-18-14-4-NESTER	.018	14	4	11.1
	G23522	MWCE-18-14-5-NESTER	.018	14	5	8.9
	G26989	MWCE-18-14-6-NESTER	.018	14	6	7.4
	G26990	MWCE-18-14-8-NESTER	.018	14	8	5.6
	G26986	MWCE-18-14-10-NESTER	.018	14	10	4.5

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.

## NOW AVAILABLE IN A BROADER SELECTION OF MACROCOILS

Rows shaded in blue are new additions to the product line.

	Order Number	Reference Part Number	Recommended Catheter ID/ End Hole Diameter inch	Extended Embolus Length cm	Coiled Embolus Diameter mm	Approximate Number of Loops
<b>.035 Inch Coils</b>						
7cm Extended Embolus	G23523	MWCE-35-7-3-NESTER	.035	7	3	7.4
	G52742	MWCE-35-7-4-NESTER	.035	7	4	5.6
	G23524	MWCE-35-7-5-NESTER	.035	7	5	4.5
	G52743	MWCE-35-7-6-NESTER	.035	7	6	3.7
	G52744	MWCE-35-7-8-NESTER	.035	7	8	2.8
	G52745	MWCE-35-7-10-NESTER	.035	7	10	2.2
	G52746	MWCE-35-7-12-NESTER	.035	7	12	1.9
	G52747	MWCE-35-7-14-NESTER	.035	7	14	1.6
	G52748	MWCE-35-7-16-NESTER	.035	7	16	1.4
	G52749	MWCE-35-7-18-NESTER	.035	7	18	1.2
G52750	MWCE-35-7-20-NESTER	.035	7	20	1.1	
14cm Extended Embolus	G23525	MWCE-35-14-3-NESTER	.035	14	3	14.9
	G26993	MWCE-35-14-4-NESTER	.035	14	4	11.1
	G23526	MWCE-35-14-5-NESTER	.035	14	5	8.9
	G26994	MWCE-35-14-6-NESTER	.035	14	6	7.4
	G26995	MWCE-35-14-8-NESTER	.035	14	8	5.6
	G26991	MWCE-35-14-10-NESTER	.035	14	10	4.5
	G26992	MWCE-35-14-12-NESTER	.035	14	12	3.7
	G52751	MWCE-35-14-14-NESTER	.035	14	14	3.2
	G52752	MWCE-35-14-16-NESTER	.035	14	16	2.8
	G52753	MWCE-35-14-18-NESTER	.035	14	18	2.5
G52754	MWCE-35-14-20-NESTER	.035	14	20	2.2	
20cm Extended Embolus	G23527	MWCE-35-20-5-NESTER	.035	20	5	12.7
	G23528	MWCE-35-20-6-NESTER	.035	20	6	10.6
	G23529	MWCE-35-20-8-NESTER	.035	20	8	8
	G23530	MWCE-35-20-10-NESTER	.035	20	10	6.4
	G23531	MWCE-35-20-12-NESTER	.035	20	12	5.3
	G23532	MWCE-35-20-14-NESTER	.035	20	14	4.5
	G23533	MWCE-35-20-16-NESTER	.035	20	16	4.0
	G52755	MWCE-35-20-18-NESTER	.035	20	18	3.5
G52756	MWCE-35-20-20-NESTER	.035	20	20	3.2	

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.

### Customer Service

**EMEA:** EDI - [www.cookmedical.com/edi.do](http://www.cookmedical.com/edi.do)  
Distributors: +353 61239240, [ssc.distributors@cookmedical.com](mailto:ssc.distributors@cookmedical.com)  
Austria: +43 179567121, [oe.orders@cookmedical.com](mailto:oe.orders@cookmedical.com)  
Belgium: +32 27001633, [be.orders@cookmedical.com](mailto:be.orders@cookmedical.com)  
Denmark: +45 38487607, [da.orders@cookmedical.com](mailto:da.orders@cookmedical.com)  
Finland: +358 972519996, [fi.orders@cookmedical.com](mailto:fi.orders@cookmedical.com)  
France: +33 171230269, [fr.orders@cookmedical.com](mailto:fr.orders@cookmedical.com)  
Germany: +49 6950072804, [de.orders@cookmedical.com](mailto:de.orders@cookmedical.com)  
Hungary: +36 17779199, [hu.orders@cookmedical.com](mailto:hu.orders@cookmedical.com)  
Ireland: +353 61239252, [ie.orders@cookmedical.com](mailto:ie.orders@cookmedical.com)  
Italy: +39 0269682853, [it.orders@cookmedical.com](mailto:it.orders@cookmedical.com)  
Netherlands: +31 202013367, [nl.orders@cookmedical.com](mailto:nl.orders@cookmedical.com)  
Norway: +47 23162968, [no.orders@cookmedical.com](mailto:no.orders@cookmedical.com)  
Spain: +34 912702691, [es.orders@cookmedical.com](mailto:es.orders@cookmedical.com)  
Sweden: +46 858769468, [se.orders@cookmedical.com](mailto:se.orders@cookmedical.com)  
Switzerland - French: +41 448009609, [fr.orders@cookmedical.com](mailto:fr.orders@cookmedical.com)  
Switzerland - Italian: +41 448009609, [it.orders@cookmedical.com](mailto:it.orders@cookmedical.com)  
Switzerland - German: +41 448009609, [de.orders@cookmedical.com](mailto:de.orders@cookmedical.com)  
United Kingdom: +44 2073654183, [uk.orders@cookmedical.com](mailto:uk.orders@cookmedical.com)

[www.cookmedical.com](http://www.cookmedical.com)

**Americas:** EDI - [www.cookmedical.com/edi.do](http://www.cookmedical.com/edi.do)  
Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335  
E-mail: [orders@cookmedical.com](mailto:orders@cookmedical.com)

**Australia:**  
Phone: +61 738411188, 1800777222, Fax: +61 738411288, 1800077283  
E-mail: [cau.custserv@cookmedical.com](mailto:cau.custserv@cookmedical.com)

