




**AVATAR<sup>®</sup>**


**PTFE Vascular Graft**

**INTERNATIONAL PRODUCT CATALOG**





AVATAR® PTFE Vascular Grafts are made entirely from PTFE. The AVATAR® graft does not have an external wrap but instead has an open pore structure averaging 20 microns internodal distance, encouraging tissue ingrowth and making it soft and easy to suture. Helix grafts incorporate an external PTFE spiral support for use in applications where kink or crush resistance is required. Tapered and step tapered grafts can allow for better management of blood flow. All AVATAR® PTFE Vascular Grafts come sterile in dual barrier tray-within-tray packaging.



(20 microns average Internodal distance)

TYPE	CODE
Center Helix	CH
Regular	RG
End Helix	EH
Full Helix	FH
Taper	TA
Taper Helix	TH
Step Taper	ST
Step Taper Center Helix	SC



(Dual barrier packaging tray-within-tray)

## PTFE Vascular Graft Product Range. (Regular Wall Thickness)

### Arterio-Venous Grafts – Standard



NO	LENGTH	DIAMETER	PRODUCT CODE
1	20cm	5mm	PVGRG0520
2		6mm	PVGRG0620
3	45cm	5mm	PVGRG0545
4		6mm	PVGRG0645



## Arterio-Venous Grafts – Center Helix (non-removable)

NO	LENGTH	DIAMETER	PRODUCT CODE
5	45cm	5mm	PVGCH0545
6		6mm	PVGCH0645



## Arterio-Venous Grafts – End Helix (non-removable)

NO	LENGTH	DIAMETER	PRODUCT CODE
7	45cm	5mm	PVGEH0545
8		6mm	PVGEH0645



## Arterio-Venous Grafts - Step-Tapered

NO	LENGTH	DIAMETER	PRODUCT CODE
9	45cm	4-6mm	PVGST4645
10		4-7mm	PVGST4745



## Arterio-Venous Grafts - Step-Tapered Center Helix (non-removable)

NO	LENGTH	DIAMETER	PRODUCT CODE
11	45cm	4-6mm	PVGSC4645
12		4-7mm	PVGSC4745



## Peripheral Vascular Grafts – Standard

NO	LENGTH	DIAMETER	PRODUCT CODE
13	20cm	4mm	PVGRG0420
14		7mm	PVGRG0720
15		8mm	PVGRG0820
16	45cm	4mm	PVGRG0445
17		7mm	PVGRG0745
18		8mm	PVGRG0845
19	60cm	5mm	PVGRG0560
20		6mm	PVGRG0660
21		7mm	PVGRG0760
22		8mm	PVGRG0860
23	75cm	5mm	PVGRG0575
24		6mm	PVGRG0675
25		7mm	PVGRG0775
26		8mm	PVGRG0875



Biomedical Innovations



**Contact:**

**BIOVIC SDN BHD**

**12A, JALAN PJU 3/40,**

**SUNWAY DAMANSARA TECHNOLOGY PARK,**

**47810 PETALING JAYA,**

**SELANGOR, MALAYSIA.**

**enquiry@biovic.com.my**

**+603 7885 0850**

**www.biovic.com.my**